Poster/abstract #346


2025 AACR Annual Meeting

A Croset¹, A Saidi², F Malvezzi¹, T Stallons², M Mettier¹, A Wong², J Rantanen¹, S Wullschleger¹, E Tselempi¹, Y Kaufmann¹, T Lekishvili¹, S Riesenberg¹, N Pina¹, J Blunschi¹, L Abduli¹, C Reichen¹, C Lizak¹, A Schatzmann², A Goubier¹, J Torgue², D Steiner¹ 1. Molecular Partners AG, Zurich-Schlieren, Switzerland; 2. Orano Med LLC, Plano, TX, USA

Introduction

- DLL3 is a promising target for radioligand therapy as it is highly upregulated in SCLC and other high-grade neuroendocrine tumors, while not expressed in healthy tissue.
- · Designed Ankyrin Repeat Proteins (DARPins) are binding proteins with high specificity and affinity that can be generated against a broad range of tumor targets and thus serve as ideal vectors.
- · By leveraging the intrinsic properties of DARPins and the learnings from our platform optimization we achieved efficient tumor uptake and penetration, while limiting exposure of healthy tissues.
- ²¹²Pb is a radioisotope with a short decay half-life and a favorable decay chain, allowing high energy deposition on tumor in a short time frame.
- Here we present preclinical results of MP0712, our DLL3-targeting ²¹²Pb-based Radio-DARPin Therapeutic (RDT) combining the advantages of a small protein-based delivery vector and the shortlived alpha particle-emitting radioisotope ²¹²Pb.

MP0712 characteristics and target indication



Small protein-based delivery vector based on Radio-DARPin platform (HLE, half-life extension)

MP0712 properties

- Specific binding with high affinity
- Affinity to hDLL3: 0.2 nM by SPR
- Human cell binding: ~ 2 nM on NCI-H82 ± HSA
- Good developability

²¹²Pb advantages

- Safety: alpha precursor with clean decay chain
- Efficacy: high energy deposition on tumor in short time (half-life of 10.6 h)


SCLC as indication

- Aggressive cancer with high unmet medical need 2L: mPFS ~3 mos; 5-year OS ~3%^{1,2}
- DLL3 is expressed in >85% of patients³

DLL3 a promising target

- Homogeneous & low tumor expression
- No expression in healthy tissues
- New treatments with room for improvement: Tarlatamab for 2L+: ORR ~40% / DoR 9.7 mos⁴

MP0712 showed good efficacy & tumor reduction



In vivo efficacy study of MP0712 in NCI-H82 tumor bearing mice


MP0712 was injected 4 x 10uCi every week or every 2 weeks (40 μCi cumulative) and 8 x 5 μCi every week or every 2 weeks (40 μCi cumulative) into the tail vein of R2G2 mice xenografted subcutaneously with NCI-H82 cells. Animals were under observation daily and 3x per week tumors were measured by caliper. Animals were sacrificed when tumors reached 1500 to 2000mm³ or if termination criteria were met. The data are expressed as average +/- SEM on the graph above and as single curve on the graphs below. Dotted lines represent the different cycle of injections (Cx-wx).


Days since cell injection MP0712 (4x10µCi, 2 wk) in NCI-H82


unnett's multiple comparisons test

Buffer only vs. MP0712 (8 x 5 μ Ci, 1 wk)

Buffer only vs. MP0712 (8 x 5 μ Ci, 2 wks)







• ~20% of mice showed complete tumor regression with MP0712 injected every week at 8 x 5 µCi, • ~20% of mice showed complete tumor regression with MP0712 injected every 2 weeks at 4 x 10 μCi

Complete tumor regression is stable after treatment MP0712 treatment resulted in complete Adjusted P-value tumor regression


Buffer only vs. MP0712 (4 x 10 μCi, 1 wk) < 0.0001 Significant effect for the doses of 4 x 10 µCi injected every week and every 2 weeks Buffer only vs. MP0712 (4 x 10 µCi, 2 wks) 0.0006

0.0002


>0.9999


Significant effect for the doses of 8 x 5 µCi injected every week but not significant when injected every 2 weeks


MP0712 showed a favorable safety profile

























- Complete recovery of body weight loss after 10 days
- Complete recovery of hematologic profile after 28 days
- MP0712 treatment up to 30 μCi well tolerated

Conclusion


- MP0712, the first ²¹²Pb-DLL3 Targeted Radio-DARPin Therapy
 - High tumor uptake
 - Reached T:K > 2 in mouse models expressing DLL3
 - Induced good efficacy & tumor reduction
 - Showed a favorable safety profile in vivo up to 40 μCi
- IND-enabling package completed
- Initial first-in-human clinical data expected in 2025


MP0712 showed favorable biodistribution and tumor specificity


















- different mouse models expressing DLL3 ²¹²Pb-DOTAM-DARPin injected 1x10µCi at 0.01mg/kg
- MP0712 reached T:K ratios >2 in mouse models matching clinically relevant
- expressing tumors confirmed high target specificity of MP0712



- - DLL3 expression levels Selective uptake in DLL3-

Acknowledgments: Quy Nguyen, Bradley Himes, Varaidzo Makina, Emily Brocious, Alyssa De Sisto, Caroline Crisp, Daniela Tapia Hernandez, and Abigail Fondren Conflict of interest disclosures: All authors with affiliation 1) are employees and stock owners of Molecular Partners; and all authors with affiliation 2) are employees of Orano Med.